

# Methods in Phase Quantification with PXRD


---

Nate Henderson, Ph.D.  
Senior Applications Scientist - XRD

# Reflection Geometry (Bragg-Brentano)




## Line Geometry (1D Detector)



- Diverging beam geometry increases sampling statistics through large sample illumination
- Rapid data collection with linear detectors (1D)
- Generally preferred over spot beam analysis (2D) for quantification


# Quantification Overview

- Concentration of crystalline phases calculated as a function of **intensity**
- Data quality is generally a limiting factor
- Selected Methods:
  - Calibration Curve
  - Reference Intensity Ratio
  - Full Pattern Summation
  - Quantitative Rietveld Refinement
- Degree of crystallinity / Percent Amorphous





# Calibration Curve

Calibration standards provide an increase in instrument response with increased concentration



# Calibration Curve

Non-reproducible preferred orientation leads to poor quality calibration curve



# Calibration Curve

---

## Benefits

- Sensitive to low concentrations (micrograms under controlled conditions)
- Simple, straight-forward concept


## Challenges

- Requires availability of standards
- Requires intensity monitoring for tube aging
- Susceptible to preferred orientation

## Reference Intensity Ratio (RIR Method)

$I/I_c$  is measured or calculated with a 50/50 wt% mixture of a given material with corundum ( $\text{Al}_2\text{O}_3$ )


$$\frac{I_a}{I_b} = \frac{(I/I_c)_a}{(I/I_c)_b} \frac{x_a}{x_b}$$





# Reference Intensity Ratio (RIR Method)

Scaling of largest reflections (100% peak) for each phase provides concentration





# Reference Intensity Ratio (RIR Method)

---

## Benefits


- Quick and simple process
- Widely implemented in XRD software packages
- I/I<sub>c</sub> values widely available or easily obtainable by user

## Challenges


- Susceptible to errors in data quality from preferred orientation and graininess
- Relies on availability of DB values for I/I<sub>c</sub> or ability of user to derive independently
- Errors with crystallite size broadening (changes to observed peak maximum)

# Full Pattern Summation

Summation of user-collected reference patterns (both crystalline and amorphous possible)



# Full Pattern Summation



# Full Pattern Summation

---

## Benefits


- Extends quantification to amorphous phases and textured phases that are highly reproducible
- Can be implemented for multiple amorphous phases

## Challenges

- Requires characterization of pure phases for baseline “fingerprint”
- Requires similar data collection parameters
- Susceptible to variations in non-reproducible preferred orientation

# Quantitative Rietveld Refinement


Standardless, least-squares modeling using scaling of calculated, known crystal structures



Hill and Howard, 1987

# Quantitative Rietveld Refinement


Refinement model can be modified and extended by adding additional structures as needed





# Quantitative Rietveld Refinement

Complex mixtures can be refined with high quality data and thorough phase identification



Phase			
	Actual (wt%)	Refined (wt%)	Difference
Quartz	15.3	16.2	0.9
K-Feldspar	8.6	8.1	0.5
Plagioclase	10.9	10.0	0.9
Calcite	10.6	12.5	1.9
Pyrite	3.9	3.6	0.3
Rutile	3.0	3.2	0.2
Anatase	1.3	1.1	0.2
Amphibole	0.1	-	0.1
Chloritoid	0.3	-	0.3
Zircon	0.1	-	0.1
Kaolinite	21.6	17.7	3.9
Muscovite	4.9	8.2	3.3
Montmorillonite	9.8	10.5	0.7
Pyrophyllite	9.6	8.9	0.7

Kaduk et al. Nature Reviews Methods Primer 2021



# Quantitative Rietveld Refinement

---


## Benefits

- Generally regarded as gold standard for quantification via PXRD
- No standards required
- Robust quantification method with ability to deconvolve overlapping reflections in complex mixtures
- Modeling possible for almost any definable parameter (e.g., preferred orientation, crystallite size broadening, anisotropic peak broadening, etc.)

## Challenges

- Requires known crystal structures for building refinement model
- Potential error due to microabsorption (contrast in densities)
- Steeper learning curve with more advanced understanding on diffraction experiment

# Degree of Crystallinity (Ratio of Integrated Areas)



# Degree of Crystallinity (Ratio of Integrated Areas)

---

## Benefits

- Simple, rapid calculation
- Widely implemented in XRD software packages
- Can analyze multiple amorphous phases in a single diffraction specimen


## Challenges

- Susceptible to user bias in integration and background fitting
- Inaccuracies possible with aggressive background modeling
- Decreased accuracies with large difference in scattering potential between amorphous and crystalline phases (different chemical compositions)

## Supplemental Reading

---

- Madsen, I.C.; Scarlett, N.V.Y.; Kleeberg, R.; Knorr, K.  
*Chapter 3.9: Quantitative Phase Analysis*  
in *International Tables for Crystallography, Volume H, Powder Diffraction*  
Eds. Gilmore, Kaduk, and Schenk, Wiley, **2019**
- Madsen, I.C. and Scarlett, N.V.Y.  
*Chapter 11: Quantitative Phase Analysis*  
in *Powder Diffraction: Theory and Practice*  
Eds. Dinnebier and Billinge, Royal Society of Chemistry, **2008**.
- Madsen, I.C.; Scarlett, N.V.Y.; Kern, A.  
*Description and survey of methodologies for the determination of amorphous content via X-ray powder diffraction.*  
Z. Krist. 226, **2011**, 944.
- Madsen, I.C.; Scarlett, N.V.Y.; Riley, D.P.; Raven, M.D.  
*Quantitative Phase Analysis using the Rietveld Method*  
in *Modern Diffraction Methods*  
Eds. Mittemeijer and Welzel, Wiley, **2012**.



Innovation with Integrity